

NEW PRODUCT HI-LIGHT

COATED ALUMINUM ORTHOPEDIC HAND & HAND FIXATION DEVICE (see reverse side)



PM-962, Large



PM-963, X-Large

Padgett Instruments® by Miltex® now includes Coated Aluminum Orthopedic Hands and Hand Fixation Device (see reverse side) for Hospital, Surgery Center and Office-based surgical procedures.

The new Coated Aluminum Orthopedic Hands are ideal for reconstructive hand surgery and provide an environmentally safe alternative to traditional lead hands and splints. Orthopedic Hands are made of pliable aluminum with a non-latex coating and are reusable and easy to sterilize. In addition, Aluminum Hands may be used during x-ray and fluoroscopy procedures.

REF	Description	Length	Width
PM-962	Coated Aluminum Orthopedic Hand, Reusable, Large	9-1/2" (24.1 cm)	8-1/2" (21.6 cm)
PM-963	Coated Aluminum Orthopedic Hand, Reusable, X-Large	13" (33 cm)	10-1/2" (26.7 cm)

See Reverse Side For Additional New Padgett Instruments®

Illustrations and content provide general description only and may be subject to change. For additional information, an in-office demonstration, or to place an order, please contact your local Miltex authorized Distributor or Miltex Customer Service at 1-866-854-8300, e-mail customerservice@miltex.com, or visit us at www.miltex.com.

Miltex and Miltex Redefining Excellence are registered trademarks of Miltex, Inc., a subsidiary of Integra LifeSciences Corporation. Padgett Instruments is a registered trademark and the Integra wave logo is a trademark of Integra Life Sciences Corporation. Copyright © 2009 Integra LifeSciences Corporation

Miltex
Redefining Excellence
INTEGRA

589 Davies Drive
York, PA 17402
phone 717 840-9335
toll free 866 854-8300
fax 866 854-8400
www.miltex.com

HAND FIXATION DEVICE



PM-961

The Hand Fixation Device is ideal for both right and left hand procedures. The Hand Fixation Device is manufactured from the highest quality medical grade German Stainless Steel and features PEEK material on the finger restraints.

- Fixation of single joints from carpus to fingertips
- No slipping or escaping of single joints
- Adjustable finger and wrist restraints make the hand fixation device convenient for use on a variety of patients
- V2A antimagnetic high-grade steel
- Fully autoclavable

REF	Description	Length	Width
PM-961	Hand Fixation Device	11-1/2" (29.2 cm)	9-1/8" (23.2 cm)