

MILTEX HI-LIGHT

BLADE-SAFE


IDEAL FOR THE SAFE REMOVAL OF
SCALPEL BLADES FROM KNIFE HANDLES



FEATURES

- Protects hands against accidental cuts and punctures
- Economically priced
- Suitable for all types of medical practices
- Easy to use
- Stainless steel
- Autoclavable

1 - Insert blade side up and align to guide



2 - Press down

3 - Pull off handle

Miltex
Redefining Excellence

589 Davies Drive
York, PA 17402
phone 717 840-9335
toll free 866 854-8300
fax 866 854-8400
www.miltex.com

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email customerservice@miltex.com, or visit us at www.miltex.com.

0103/2/B/1

SURG-BLD-SAFE


H834SURGBLDSAFES