


MILTEX HI-LIGHT

FIBER OPTIC LIGHT SOURCE

150-watt Xenon



Miltex offers a complete line of state-of-the-art illumination products and accessories for the operating room and for office-based surgical procedures.

The 150-watt Xenon fiber optic light system provides superior illumination for optimal visualization. This system features a multiport turret for easy selection of various cable connection; including, Olympus, Storz®, Wolf and ACMI-styles. An adjustable attenuator knob enables the user to select the desired light intensity from 0 - 100%. Additional features include a cumulative lamp hour usage meter and easy access for lamp replacement which does not require special tools.

The Xenon lamp is environmentally friendly and unlike metal halide lamps, the Xenon lamp does not require special waste handling for mercury.

See reverse side for additional information on accessories such as Floor Stand, Autoclavable Fiber Optic Cable, Retractors, etc.

REF#	DESCRIPTION
EM-2151LS	Fiber Optic Light Source, 150-watt Xenon

SPECIFICATIONS

DESCRIPTION:	150-watt Xenon light source with selectable multiport wheel (Standard Ports: Olympus, Storz, Wolf and ACMI-styles)
CONTROLS AND INDICATORS:	Main and Standby Power Switches, Cumulative Lamp Age Indicator, Lamp On Indicator, Lamp Brightness Control 0-100%
DIMENSIONS:	14.0"L x 9.3"W x 4.3"H (356mm L x 236mm W x 110mm H)
WEIGHT:	12.2 lbs. (5.53kg)
POWER INPUT:	100-120 VAC, 200-240 VAC 50-60 Hz ± 10%
USER MANUAL:	REF# IFU-EM2151LS

See reverse side for additional information

Illustrations and content signify general description only and may be subject to minor changes



+H834SURGFOLIGHTHI05


Miltex
Redefining Excellence

589 Davies Drive
York, PA 17402
phone 717 840-9335
toll free 866 854-8300
fax 866 854-8400
www.miltex.com

ACCESSORIES



EM-2151STAND



EM-2151BRACK



EM-2151



EM-2155

ACCESSORIES

REF#	DESCRIPTION
EM-2151STAND	Fiber Optic Light Source Floor Stand w/Handle for EM-2151LS
EM-2151BRACK	Fiber Optic Light Source Bracket Adapter for EM-2151STAND
EM-2151	Fiber Optic Light Cable, 8' (244 cm), Autoclavable, with universal Olympus & ACMI male light source connection and female ACMI Instrument Coupler
EM-2155	Fiber Optic Cable Adapter for Storz and male ACMI



PM-884FO



PM-529FO



PM-530FO



PM-895FO



PM-890FO



PM-5551FO



PM-490FO
PM-491FO



PM-4570FO
PM-4571FO



PM-9057FO



PM-4270FO

FIBER OPTIC RETRACTORS (All come with standard fittings including, Storz, Olympus, Wolf and ACMI)

REF#		Handle Length	Blade Width
PM-884FO	Biggs Face Lift Retractor	9" (229 mm)	5.5 cm
PM-529FO	Biggs Mammoplasty Retractor, Standard	8" (203 mm)	5 cm
PM-530FO	Biggs Mammoplasty Retractor, Narrow	6" (152 mm)	3 cm
PM-895FO	Ferreira Breast Augmentation Retractor	5 1/8" (130 mm)	1" (25 mm) x 6 3/4" (171 mm)
PM-890FO	Ferreira Face Lift Retractor	4-1/4" (108mm)	7/8" (22 mm) x 4 1/4" (108 mm)
PM-5551FO	Aufricht Retractor	8" (203 mm)	7 mm
PM-490FO	Freeman Flap Retractor	8" (203 mm)	1/2" (13 mm)
PM-491FO	Freeman Flap Retractor	8" (203 mm)	1" (25 mm)
PM-4570FO	Deaver Retractor	8" (203 mm)	1/2" (13 mm)
PM-4571FO	Deaver Retractor	8" (203 mm)	1" (25 mm)
PM-9057FO	Padgett-Deaver Retractor	12" (305 mm)	1" (25 mm)
PM-4270FO	Killian Speculum	5 1/2" (140 mm)	2" (51 mm)

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex®, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email customerservice@miltex.com, or visit us at www.miltex.com.