


MILTEX® HI-LIGHT

ATRAUMATIC TENACULUM FORCEPS SCHROEDER & BRAUN



30-966ATR


Standard



Schroeder

Braun

New



Schroeder
30-966ATR

Braun
30-965ATR

Traditional patterns of Tenaculum Forceps have overlapping tips which can cause trauma and bleeding when used to firmly grip the Cervix, Uterus and other general tissue applications. The new Miltex Atraumatic Tenaculum Forceps are designed with non-overlapping tips, even when closed with the ring handle ratchet fully engaged. The tip design insures a firm and reliable grip of tissue, but minimizes patient bleeding and discomfort. A "must have" instrument to complement your procedural sets for Colposcopy, Endometrial Biopsy and other commonly performed Gynecological procedures.

The Miltex name is synonymous with premium surgical instruments, manufactured from the highest quality stainless steel German forgings and crafted to exacting specifications.

FEATURES & BENEFITS:

- *Atraumatic design with non-overlapping tips even when fully closed.*
- *Allows for gentle grasp and hold of tissue during cervical, uterine and other pelvic procedures.*
- *Ring-handles are finished with rounded edges for optimum grip comfort.*
- *Ratchet design provides smooth and controlled action*
- *Jaw options include either: Rounded (Schroeder) and Straight/Squared (Braun)*
- *Available in both Miltex and MeisterHand*

REF	BRAND	DESCRIPTION
30-965ATR	Miltex®	Braun Tenaculum Forceps, 10" (25.4 cm), straight/squared jaw, with non-overlapping atraumatic tips
30-966ATR	Miltex®	Schroeder Tenaculum Forceps, 10" (25.4 cm), rounded jaw, with non-overlapping atraumatic tips
MH30-965ATR	MeisterHand™	Braun Tenaculum Forceps, 10" (25.4 cm), straight/squared jaw, with non-overlapping atraumatic tips
MH30-966ATR	MeisterHand™	Schroeder Tenaculum Forceps, 10" (25.4 cm), rounded jaw, with non-overlapping atraumatic tips

In addition to our newly featured Atraumatic Tenaculum Forceps, Miltex offers a comprehensive offering of other OB/GYN instruments and disposable products to meet a variety of your clinical needs.

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email customerservice@miltex.com, or visit us at www.miltex.com.



+H834SURGTENACHIO\$

Miltex
Redefining Excellence

589 Davies Drive
York, PA 17402
phone 717 840-9335
toll free 866 854-8300
fax 866 854-8400
www.miltex.com