


MILTEX® HI-LIGHT

NEW! SIMS AND SIMPSON UTERINE SOUNDS MALLEABLE STAINLESS STEEL



Miltex now offers German-crafted SIMS and SIMPSON Uterine Sounds manufactured from malleable and fully autoclavable stainless steel. The new bendable stainless steel uterine sounds offer clinicians a reliable, quality instrument for both routine office-based procedures or in the operating room. The unique stainless steel design is a helpful alternative to traditional patterns manufactured with silver or chrome plating which may often contribute to tarnishing and/or stain deposit transfer during autoclaving and cleaning processes. The new Miltex flexible stainless steel uterine dilators reduce time separating and processing dissimilar metal instruments and eliminates scrubbing off instrument stains.

REF	Description	Length
30-650SS	SIMS Uterine Sound	13" (33 cm), 4mm OD, graduated in cm, malleable, Stainless Steel
30-655SS	SIMPSON Uterine Sound	13" (33 cm), 4mm OD, graduated in inches, malleable, Stainless Steel

Illustrations and content signify general description only and may be subject to minor changes. For additional information, an in-office demonstration, or to place an order, please contact your local Miltex authorized Distributor or Miltex Customer Service at 1-866-854-8300, e-mail customerservice@miltex.com, or visit us at www.miltex.com.

Miltex
Redefining Excellence

589 Davies Drive
York, PA 17402
phone 717 840-9335
toll free 866 854-8300
fax 866 854-8400
www.miltex.com

