


# MILTEX® HI-LIGHT

## NEW AND IMPROVED McGIVNEY HEMORRHOIDAL LIGATOR

(FOR USE WITH BOTH LATEX-FREE & STANDARD LATEX O-RINGS)



28-154B (includes loading cone REF# 28-156)

28-155

28-152

### PRODUCT DESCRIPTION

Improved lightweight model with comfortable offset handle for controlled grip and better viewing access. The shaft with ligating drum can be rotated 360 degrees for correct positioning and is easy to disassemble for thorough cleaning. Stainless Steel loading cone facilitates placement of the O-ring bands onto the ligating drum. Compatible with Latex and Latex-Free ligating O-rings. (See reverse side for additional product details.)

#### REF

<b>28-154A</b>	McGIVNEY Hemorrhoidal Ligator Set for LATEX-FREE applications - includes McGivney Hemorrhoidal Ligator 8" shaft (PN 28-154B), McGIVNEY Grasping Forceps 7-1/2" (PN 28-152), Loading Cone (PN 28-156) and 100 LATEX-FREE O-Rings (28-155).
<b>28-154B</b>	McGIVNEY Hemorrhoidal Ligator improved model for LATEX-FREE applications - 8" (20.3 cm) working length. Supplied with Loading Cone (28-156) and 100 LATEX-FREE O-Rings (28-155).
<b>28-152</b>	McGIVNEY Hemorrhoidal Grasping Forceps, 7-1/2" (19.1 cm), angled shanks. Specially designed for use with McGivney Ligators.
<b>28-155</b>	LATEX-FREE O-Rings compatible with PN 28-154B Ligator only. 100 rings per bag. Provided Non-Sterile.
<b>28-151</b>	LATEX O-Rings compatible with PN's 28-150, 28-150XL & 28-154B Ligator patterns. 100 rings per bag. Provided Non-Sterile. (See reverse side for additional product details.)
<b>28-156</b>	Loading Cone compatible with PN 28-154B Ligator only.

### See reverse side for additional product details

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email [customerservice@miltex.com](mailto:customerservice@miltex.com), or visit us at [www.miltex.com](http://www.miltex.com).




589 Davies Drive  
York, PA 17402  
phone 717 840-9335  
toll free 866 854-8300  
fax 866 854-8400  
[www.miltex.com](http://www.miltex.com)



# McGIVNEY HEMORRHOIDAL LIGATOR

(FOR USE WITH LATEX O-RINGS ONLY)



28-150 (includes loading cone REF#: 28-153)

28-151

28-152

## PRODUCT DESCRIPTION

Improved lightweight model with comfortable offset handle for controlled grip and better viewing access. The shaft with ligating drum can be rotated 360 degrees for correct positioning and is easy to disassemble for thorough cleaning. Stainless Steel loading cone facilitates placement of the O-ring bands onto the ligating drum. Compatible with Latex ligating O-rings only.

### REF

<b>28-150</b>	McGIVNEY Hemorrhoidal Ligator, 7" (17.8 cm) working length. Supplied with Loading Cone (28-153) and 100 LATEX O-Rings (28-151)
<b>28-150XL</b>	McGIVNEY Hemorrhoidal Ligator Extra-Long, 10" (25.4 cm) working length. Supplied with Loading Cone (28-153) and 100 LATEX O-Rings (28-151)
<b>28-151</b>	LATEX O-Rings compatible with PN's 28-150, 28-150XL & 28-154B Ligator patterns. 100 rings per bag. Provided Non-Sterile.
<b>28-152</b>	McGIVNEY Hemorrhoidal Grasping Forceps, 7-1/2" (19.1 cm), angled shanks. Specially designed for use with McGivney Ligators.
<b>28-153</b>	Loading Cone compatible with PN 28-150 & 28-150XL Ligators only.